

## Supplement kit for microvascular endothelial cell growth medium, advanced, HuS-kit

Cat.-Nr.: 216 1102

contains of:

Basal media	Supplements
-	222 1000 L-Glutamine 232 2500 HuS (human serum AB) 244 0250 human rec. EGF (epidermal growth factor) 245 0500 human rec. bFGF (basic fibroblast growth factor) 242 1000 human rec. Long R3 IGF-1 223 0005 Ascorbic acid 241 0025 human rec. VEGF 224 0050 Hydrocortisone 236 0350 Antibiotics (optional)

### Maintenance of supplement kit for microvascular endothelial cell growth medium, advanced:

The supplements are delivered on dry ice. Store the **supplements** at **-20°C**.

### Characteristics:

The Provitro supplement kit for microvascular endothelial cell growth medium, advanced is suitable for culturing Provitro human microvascular endothelial cells **after adding the supplement kit components to 500 ml basal medium**. The formulation is optimized for initial seeding of 4,000 cells / cm<sup>2</sup> up to confluence (approx. 90 %). Feeder-layer, matrix substrates or other substances are not necessary. **Due to the possibility of reduced proliferative activity we recommend to use the antibiotic supplement for freshly isolated cells only.**

### Stability and storage:

Store the **supplements** at **-20°C**. Do not use the supplements after expiring date.

### Special note:

Avoid any thaw/freeze cycle.

### Quality control:

Provitro's supplement kit for microvascular endothelial cell growth medium, advanced is thoroughly tested after each production. All components are tested in a stringent biological assay. Each batch is checked for HMVEC proliferating characteristics. The cells cultured in microvascular endothelial cell growth medium are checked regarding their morphology, the adherence rate, the colony forming efficiency and the population doubling time.

### In vitro laboratory use only.

Not intended for any human or animal diagnostic or therapeutic use.